

UNDERWATER BICYCLE

NOT YOUR ORDINARY BICYCLE

BUT STRIKINGLY SIMILAR



EWAC
MEDICAL
We get you moving

www.ewacmedical.com



We believe that in rehabilitation, you should be able to offer a wide range of exercise options. We also believe that exercising in water should be fun. The EWAC Underwater Bicycle is designed for cycling exercises in water for rehabilitation purposes. It is equipped with a unique resistance mechanism, which causes resistance to increase exponentially with the cycling speed. The frame is made of polished 316 stainless steel and therefore suitable for use in a therapy pool.

sensation the bike is still going. The resistance increases exponentially with velocity.

Exercising in water increases the efficiency of the heart and because of the heat conductivity, the body is easily cooled while the connective tissue becomes more flexible. Water is a safe and low impact exercise environment. The underwater bicycle is an excellent tool for rehabilitation as well as combating obesity.

The unique resistance mechanism closely resembles the conditions actual cycling. Increasing speed takes a considerable effort. Slowing down, the inertia gives the

EWAC Medical is Lloyd's Register Quality certified according to ISO 9001 and ISO 13485.



EWAC Medical
Galileistraat 20
1704 SE Heerhugowaard
The Netherlands

**EWAC
MEDICAL**
We get you moving

+31 226 313 457
medical@ewac.nl
www.ewacmedical.com

Subject to alterations and errors. No rights may be derived from this brochure