


Transfer resources


Pool chair

- seat and back HT upholstery
- hoisting eyes on back rest and front legs
- separately foldable armrests, ergonomically shaped
- armrests in aluminium, white coated
- four stainless steel swivel castors, Ø 125 mm with double brake
- adjustable and foldable footrests
- only available in seat height 490 mm
- only available in stainless steel 316 polished version
- customization not possible



Swimming pool chair self-propelled

- seat and back HT upholstery
- suspension eyelets on backrest and front legs
- separately foldable armrests, ergonomically designed
- armrests in aluminium, white coated
- stainless steel swivel castors Ø 125 mm with double brake and plastic 24" castors with brake and hoop
- adjustable and foldable footrests
- only available in seat height 490 mm
- only available in stainless steel 316 polished version
- customization not possible



Pool stretcher

- lying surface HT-upholstery
- four stainless steel swivel castors, Ø 125 mm with double brake
- Adjustable backrest by means of gas spring
- height adjustable by means of hydraulic pump with left and right operation
- two folding side rails, which are mounted on the frame insert
- Upper frame is removable from lower frame and equipped with lifting eyes
- only available in the indicated sizes with the length 1900 mm
- only available in stainless steel 316 polished execution
- customization not possible
- maximum user weight 150 kg.

