## WATER PERMEABLE SURFACE

Pool water can move **through the surface** undisturbed, because of openings between the tiles.Therefore:

- Optimal water hygiene
- Undisturbed water circulation
- The floor moves faster through the water

## WATER HYDRAULICS



Water hydraulic cylinder means zero risk of contamination with toxic fluids

#### MINIMUM CONSTRUCTION THICKNESS



A total construction thickness of 30 cm, this is the thinnest movable floor in the market.

## MAXIMUM ADAPTABILITY

- Underwater treadmill
  - Underwater bicycle



- Camera's
- Counter current

Underwater massage

Self adjusting steps

## QUALIFICATIONS

- Certified as Class I Medical Device
- ISO 9001 Certified manufacturer
- ISO 13485 Qualified supplier

## **RETRO FITTING POSSIBLE**

Can be installed in new and existing pools

#### SAFE

- Strong stainless steel frame
- Safe non-slip surface
- Water permeable over entire surface
- Safe water hydraulics

## DYE TEST

Movable floor that passes the EN15288-2:2018 annex A Dye test

# Additional advantages





years

