



# MOVABLE FLOOR


## WATER PERMEABLE SURFACE

Pool water can move **through the surface** undisturbed, because of openings between the tiles. Therefore:

- ♥ Optimal water hygiene
- ➔ Undisturbed water circulation
- The floor moves faster through the water




## WATER HYDRAULICS



**Water hydraulic cylinder** means zero risk of contamination with toxic fluids

## MINIMUM CONSTRUCTION THICKNESS



A **total construction thickness of 30 cm**, this is the thinnest movable floor in the market.

## MAXIMUM ADAPTABILITY

- Underwater treadmill
- Underwater bicycle
- Obstacle course
- Camera's
- Counter current
- Underwater massage
- Self adjusting steps

## SAFE

- Strong stainless steel frame
- Safe non-slip surface
- Water permeable over entire surface
- Safe water hydraulics

## DYE TEST

Only movable floor that passes the EN15288-2:2008 annex A Dye test

## QUALIFICATIONS

- Certified as Class I Medical Device (MDD)
- ISO 9001:2008 Certified manufacturer
- ISO 13485:2012 Qualified supplier

## RETRO FITTING POSSIBLE

Can be installed in new and existing pools

## Additional advantages

Ready for use in 1 week

Optimal water hygiene

**25+** years lifespan

Single point of contact

Easy maintenance

**80%** recyclable materials

No toxic fluids

**800+** projects worldwide