











# POOLTRACK®

Underwater Treadmill



Informative and programmable  
touchscreen interface

## POOLTRACK® UNDERWATER TREADMILL

-  Perform stationary walking exercises
-  Study gait and walking pattern with stationary camera's
-  For rehabilitation, physical training or weight loss
-  Safe 17V AC drive system
-  Speeds from 0,2 to 11 km/h
-  Use the same space for other purposes when not in use
-  Intervals programmable
-  No toxic fluids
-  Wide walking surface
-  Low maintenance



EWAC Medical  
Dulleweg 18  
1721 CT Langedijk  
The Netherlands

**EWAC  
MEDICAL**  
*We get you moving*

+31 226 313 457  
medical@ewac.nl  
www.ewacmedical.com