

WATER TREATMENT



WATER TREATMENT

Apart from state of the art Aquatic Therapy equipment, EWAC Medical also takes care of the entire water treatment plant to ensure optimal hygienic safety for the patient.



WATER TREATMENT

Using modern and well tested technologies, EWAC Medical builds the right system that suits the size and use of the pool. Especially in medical treatment situations, hygienic safety of the water is paramount. With our water treatment systems the chlorine concentration can be kept at a value that is both low and safe, and also the pH is kept at a safe and comfortable level for the users.



t = 0



t = 15



DYE TEST

The dye test is described in annex A of the EN 15288-2:2008. A dye test will prove that no 'dead pockets' are present in the swimming pool. EWAC medical's movable swimming pool floor is the first in Europe which has successfully passed a dye test.

WATER TREATMENT

In combination with the EWAC Medical movable floor system, our systems are unique in the market because they conform to the strict 15288-2:2008 Annex A Dye test. This means that regardless of the presence of an adjustable floor in the water, the water is still properly circulated and bacteria are still effectively countered by the water treatment system.